

The Original **MOT Foot Mold System**




Create Lightweight Durable Positive Molds

- High Quality • Easy to Use • Precise
- Time Savings • Material Savings



CASTING & FABRICATION

Provides flat bed and/or Mandril forming capabilities common to the majority of fabricator's techniques.



Use MOT Rigid Foam #1, #2 or #3 or your normal mold medium.

MOT RIGID FOAMS are easily modified on a trautman carver or by hand. Plaster modifications are compatible. MOT RIGID FOAMS are strong, lightweight two part urethane foams that are cost effective and very simple to use.



Reduce Mold Breakage, Material & Fabrication Costs



MYRDAL
Orthopedic Technologies

For more information contact **Myrdal Orthopedic Technologies:**
TOLL FREE FAX: 866-522-2722 • TOLL FREE PHONE: 877-395-0081
Phone: 204-475-3799 • Fax: 204-775-2252
Email: mot@escape.ca • Web Site: www.myrdalorthopedics.com