

QUICK-SIL™

5 - Minute Silicone System

- Super soft, 2⁺ shore A blue silicone
- Provides a custom layer of skin compatible, viscoelastic silicone
- Very fast set time
- Self-leveling
- With Quick-Prime or Fix-Sil, will adhere to most materials
- Water permeable for dispersal of moisture
- Non-hazardous: Does not react with tissue



50ml. Starter Kit
Product #: 55/1250-005

Orthotic Work or Small Modifications

50 ml. Starter Kit Includes:

- (2) 50ml. cartridges
- Mixing Gun
- (10) static mixers
- Quick-Prime Kit
- Fix-Sil silicone adhesive

Suggested Uses:

- Pressure relief in foot orthoses
- Posterior calcaneal relief
- Ulnar Styloid relief
- Body Jacket relief



330ml. Starter Kit
Product #: 55/1250-033

Prosthetic Work and Larger Modifications

330 ml. Starter Kit Includes:

- (1) 330ml. cartridge
- Mixing Gun
- (10) static mixers
- Quick-Prime Kit
- Fix-Sil silicone adhesive

Suggested Uses:

- Distal end pads
- Forefoot prostheses and "bumpers"
- Pressure relief in seating systems
- Impact and pressure relief in sports equipment
- Myoelectric isolation



To avoid contamination, store unused Quick-Sil with static mixer attached.

UCO

INTERNATIONAL

800-541-4030

www.ucointernational.com



To slow down Quick-Sil's cure time, cool cartridges in refrigerator before use.

Quick-Sil 50ml.:

- 2" Shore A 5-pack cartridges
- Static mixers 25 Static mixers

Item	Product Number	Color	Unit
Quick-Sil 50ml. 2" Shore A Cartridges	55/1251-050	Blue	5-Pack
Quick-Sil 50ml. Static Mixers	55/1251-051	N/A	25 Pieces

Quick-Sil 330ml.:

- 2" Shore A cartridge
- Static mixers 15 Static mixers

Item	Product Number	Color	Unit
Quick-Sil 330ml. 2" Shore A Cartridges	55/1251-330	Blue	Cartridge
Quick-Sil 330ml. Static Mixers	55/1251-331	N/A	15 Pieces

Quick-Prime:

- Kit Catalyst, Base, Palette Knife, Cotton Swabs, Mixing Plate
- Components Catalyst and Base

Item	Product Number	Color	Unit
Quick-Prime Kit and Components	55/1251-440	N/A	Kit and Components
Quick-Prime Components	55/1251-444	N/A	Components

Fix-Sil Silicone Adhesive:

Fix-Sil is an effective way to repair damage to silicone products. It works to bond silicone to silicone and other materials.

Instructions:

Place small amount along tear line

Suggested Uses:

- Quick-Sil reattachment
- Repair silicone liners
- Silicone heel and insole repair

Item	Product Number	Color	Unit
Fix-Sil Silicone Adhesive	55/1251-442	Clear	3 oz. Tube



Use Fix-Sil to attach top layers to Quick-Sil modifications.

Distributed by: **Vittoria Phoenix inc**



70 Brock Street,
Vittoria, ON N0E 1W0
Tel: (800) 563-1296
Fax: (888) VP FAX IT
vphoenix@on.aibn.com
www.orthoVP.com



INTERNATIONAL
800-541-4030

