

# VPhoam

VPhoam is a cross-linked Polyethylene and EVA foam of a very high quality. It is manufactured in Europe under ISO standards.

Acceptance of this material has been overwhelming. For that reason we are now stocking three colours: Black, Blue and Flesh.

We have introduced Perforated versions in 2mm and 3mm thickness in Black, Blue and in Flesh.

Users like the fine cell structure, long lasting resilience ..... and the low price.

1.5 mm (1/16")	Black, Blue
2.0 mm (3/32")	Black, Blue, Flesh
3.0 mm (1/8")	Black, Blue, Flesh
4.0 mm (5/32")	Black
5.0 mm (3/16")	Flesh
6.0 mm (1/4")	Flesh

## Characteristics:

Density:

Shore "A" 35

Working Temperature

110 to 130 C -- 230 to 266 F

Other:

Extremely high compression resistance and very favourable compression set ratio, does not shrink, may be remolded with heat.

Suggested use:

As a cushion material for all kinds of orthoses and prostheses, soft sockets and for producing foot orthoses. An excellent alternative to Multiform™ Or Microcel Puff™.

All sheets are 1 meter x 2 meters (39" x 78")



Vittoria Phoenix inc

1534 Old Brock Street, Vittoria, ON N0E 1W0

Tel: (800) 563-1296

Fax: (888) VP FAX IT

Visit our on-line catalogue

[www.orthoVP.com](http://www.orthoVP.com)

[vphoenix@on.aibn.com](mailto:vphoenix@on.aibn.com)