


World-Class Orthotic Fabrication- Your place or ours?

In-Office Fabrication



With our VF-300™ orthotics production system, you can cast, scan, make positive casts, press the shells, and finish a pair of functional foot orthotic in less than 20-minutes per pair, right there in your own clinic. And, what's more, you can do it without plaster casts or foam impression boxes, without plaster positives and without a dirty and noisy milling machine. We're not talking pre-forms, customized insoles, or milled EVA arch

supp)™ system produces laboratory-quality (in fact, better in most cases), truly functional foot orthoses. Unlike other CAD/CAM systems that can only mill EVA, the VF-300™ uses advanced robotic-arm pin matrix technology to create a positive cast of the foot in less than 4-minutes! And because the orthotic shells are pressed not milled, you're not limited to making bulky EVA arch supports - you can press any shell material, including the advanced graphite composites.



High-end Lab Services



If you don't want to fabricate your own orthoses, you can send us plaster casts or foam boxes; or if you want to save time, mess, and cost, you can use our reusable pin casting matrix and Accu-Pin™ foot scanner (foam box re-placement), or our Accu-Light™ foot scanner (suspension plaster cast replacement) to upload "casts" to the 3-D Analytics™ software to create your own modified positive casts, and then transmit your digital images

electronically to our in-house Laboratory. We use our own VF-300™ robotic forming technology to produce accurate 3-D positive casts of your digital images, over which we can manufacture a comprehensive selection of the best foot orthoses available.



For a free demo
Call us today at 877-796-7496

www.varifitorthotics.com

7 Steps in 20 minutes!

1.



Pin Cast

1.



Light Scan

2.



Scan

3.




Educate

4.




Correct Cast

5.



Set Positive

6.



Form

7.



Finish



VF-300



Cart Scanner

 **VARIFIT™**
ORTHOTICS

Varifit Orthotics
1693 W. American Way
Suite 3
Payson, UT 84651
Toll Free Phone: (877)796-7496
www.varifitorthotics.com

Distributed by



Vittoria Phoenix
1534 Old Brock St.,
Vittoria, Ontario N0E 1W0
Toll Free Phone: (800)563-1296
Toll Free Fax: (888)873-2948
www.orthovp.com