

Feature Product



Model 328

WalkWell Universal Post Operative Shoe



WALKWELL

WalkWell's commitment to provide superior post operative footwear continues with the WalkWell Model 328.

While WalkWell has always recommended wearing matching shoes to improve balance, stability and reduce low back pain, now you can buy WalkWells in single universal units.

Of course, they wouldn't be WalkWells

without the Rigid-Rocker sole, which off-loads pressure in the metatarsal region of the foot. Or, the unique adjustable heel closure and full heel counter to fit almost any size or shape heel.

The insole, which is made of bacteria static EVA, absorbs shock and tenderly cushions the foot while the upper, made with a soft, light-weight, expandable and breathable micro-fabric, utilizes three Velcro closures to gently accommodate swelling and bandages.

They also feature a distinctly lower, slightly wider heel.

Not convinced? A study published in the Journal of the American Podiatric Medical Association indicates that WalkWells can reduce peak pressure up to 20% over

To encourage patient compliance, WalkWells were designed to be as fashionably neutral as possible - without compromising medical necessity.



The WalkWell Model 328. A universally better shoe.



WalkWells fit almost any size foot and adjust to accommodate dressings.

conventional post-operative shoes.

WalkWell Model 328s are sold in single-unit re-useable shoe bags and come in seven sizes from women's small (5) to men's xx-large (16.5).

Indications:

- Osteotomies with and without fixation
- Metatarsal fractures
- Metatarsophalangeal joint injuries
- Metatarsophalangeal tenosynovitis and tendinitis
- Trauma to the midfoot or metatarsal region
- Hallux limitus
- Hallux valgus
- Neurophathic ulceration
- Persistent edema

Vittoria Phoenix inc
70 Brock Street
Vittoria, ON N0E 1W0

Tel: (800) 563-1296

Fax: (888) VP FAX IT

www.orthoVP.com

vphoenix@on.aibn.com

Say Ahhh.

PRESENTING THE
WalkWell 502™
Adjustable Wound Care Shoe



Announcing the re-invention of the wound care shoe. The WalkWell Model 502. Universal Adjustable Wound Care Shoe, the only shoe of it's type with the acclaimed WalkWell RigidRocker sole.

The WalkWell 502 offers a wider version of the famous RigidRocker sole to accommodate patients with ulcerations or other metatarsal wounds. The RigidRocker sole offloads metatarsal pressure up to 20% more than conventional flexible wound care shoes.

Each WalkWell 502 comes with two removable, moldable and modifiable insoles, made of 1/4" Plastazote and 1/4" Poron.

This combination of materials utilizes the efficacy of Plastazote to increase the surface area and reduce pressure on the foot, contributing to the reduction of callus formation, and the unique ability of Poron to dramatically reduce vertical vector forces and plantar pressures while providing cushion and shock absorption.

The second insole enables the physician to provide patients with a fresh insert for follow-up visits.

And these WalkWells are constructed with extra depth to accommodate supplementary insoles, as necessary. Additional inserts can be purchased separately.



For flexibility, the WalkWell 502 outer shell is made with neoprene, while a soft fleece lining cushions the foot and adds absorbancy. WalkWell EasyPull Tabs simplify fitting and a padded heel and ankle closure snugly secures the foot and reduces sheer stress.